Gebrauchsanweisung\* - Informazione d'uso  
Instruction manual - Instrucciones de uso  
Mode d'emploi

Deutsch

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Date of update: 12.12.11

**Intended use:**

After being modelled and adapted, the product is inserted into the human body by specialised technicians. It is used to model dental prosthesis substituting worn out or defective anatomic parts. This product is not suited to get in contact with human tissues, bones or dentine.

**Manual for: ICE Zirkon and Prettau Zirkon**

1. Glue zirconia block into relevant insert on left side of milling table
2. Glue template with mock-up frame into relevant insert on right side of milling table.
3. Start milling process
  - a) rough milling (Milling bur 4L - 4 mm)
  - b) copy - milling (Milling bur 2L - 2 mm)
  - c) copy - milling (Milling bur 1L - 1 mm)
  - d) copy - milling (Milling bur S - 0.5 mm)
4. Remove milled frame work from block and pre-trim with diamond disc
5. Finish with tungsten cross-cut burs
6. De-dust frame with air gun
7. Colour frame with COLOUR LIQUID (see manual for COLOUR LIQUID)
8. Final-dry frame under infrared lamp for 45 minutes
9. Load furnace
10. Shut gate and select program. (see furnace manual)
11. Remove sintered framework from furnace at room temperature.
12. Check frame on working model for seat and margin fit. Trim with water cooled air turbine if required.
13. Blast using 120 micron Alu-oxide at pressure 4 bar.
14. Clean frame with brush using water. Do not steam-clean.
15. Frame is now ready for veneering

**Safety instructions:**

Ensure adequate ventilation.

For the porcelain veneering of bridges of various elements you have to observe the following parameters during the heating- and the cooling-phase:

0 - 1 g per dental-unit	55° / min.
1 - 2 g per dental-unit	45° / min.
2 - 3 g per dental unit	35° / min.
more than 3 g per dental unit	25° / min.

**Wear protective glasses.**

**Protect respiratory system:**

Wear suitable face mask to prevent ,fine dust' inhalation. Avoid ,fine dust' air pollution by using the air blaster in short bursts only rather than a continuous flow.

**Storage and Handling:**

Keep container dry. Keep at temperatures: < 28°C / < 82°F  
CAUTION: US Federal law restricts this device to sale by or on the order of a licensed dentist.

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**Emergency Phone:** +39 0474 066 660 -  
+39 348 120 15 34 (24 hours/7 days)

ICE ZIRKON and PRETTAU ZIRKON have especially been designed for the use with ICE ZIRKON MALFARBEN, ICE ZIRKON KERAMIK and COLOUR LIQUID.

**N.B. Medical products are non-returnable.**

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