

THE JAWALIGNER SYSTEM

Plaster-free articulation of dental models - Easy to design, easy to produce

Zirkon zahn



Driven by passion we have been working tirelessly over the past years to provide our customers with new solutions which ensure better results. With our technologies, we want to further improve the efficiency and quality of the daily workflows in clinics and laboratories. For this reason, we developed our new JawAligner – an innovative system to increase the techno-clinical cooperation.

furico Steque Zlian Stegni




THE JAWALIGNER

The new JawAligner are magnetic spacer plates which allow to fix in the articulator printed or milled models without using plaster. Available in different heights, they permits material saving and the articulation of models indipendently of their height. Moreover, the two sets – one inclined and one straight – allow to articulate models according to the occlusal plane inclination. Thanks to the standardised connection and the corresponding JawAligner, milled or printed models can be magnetically fixed in Zirkonzahn's PSI articulator, GSI plaster articulator and Mini-Arti ZSI to simulate patient's jaw movements. With the corresponding accessories, the JawAligner are also compatible with articulators of other manufacturers.



AN EVEN MORE COMPLETE DIGITAL WORKFLOW

With the Detection Eye intraoral scanner, the P4000 system for 3D printing and the JawAligner, we provide dentists and dental technicians with a complete pre-set package where all components are perfectly compatible with each other to ensure a simple and uniform workflow. Based on the oral situation acquired with the intraoral scanner, models are designed in the Model Maker module of the Zirkonzahn. Modifier software and printed by means of the pre-configured P4000 system for 3D printing or milled with Zirkonzahn's milling units. The models can be then plaster-free articulated into the PS1 articulator or Mini-Arti ZS1 thanks to the new JawAligner PS1 or ZS1.



Zirkonzahn

WHY USE THE JAWALIGNER?

The hardening process of plaster might cause deformations to the printed models. For this reason, a plaster-free articulation with the JawAligner is recommended in order to avoid inaccurate results. Furthermore, the use of the JawAligner ensures a precise occlusion check, which in conventional articulation methods can be compromised due to plaster expansion. In a dentistry field more and more focused on intraoral scans, the JawAligner system also allows to easily transfer the virtual patient's jaw position into the laboratory articulator for final control or case delivery and discussion.



TRADITIONAL ARTICULATION VS. JAWALIGNER SYSTEM



Long waiting times



Efficiency and flexibility: models can be designed using intraoral scans or digitised tradional impressions and then virtually articulated with the JawAligner, saving a lot of time. In addition, plaster models no longer need to be created



Full digital design



Material waste



No material wastefulness thanks to the different JawAligner heights



Possibility to create the model with a honeycomb structure ensuring greater stability



Risk of errors due to plaster expansion during the hardening phase, leading to a possible wrong occlusion



Any deformation of resin models or errors due to plaster expansion is avoided, ensuring a correct occlusion



Thanks to the two available sets, dental technicians can design models according to the patient's occlusal plane inclination, avoiding deformations caused by unbalanced front-rear height



COMPATIBILITY OF THE JAWALIGNER CONNECTION

Thanks to the standardised connection and the corresponding JawAligner, milled and printed models can be fixed in Zirkonzahn's PS1 articulator and Mini-Arti ZS1 to simulate patient's jaw movements. Optionally, it can also be used with the GS1 plaster articulator.





PS1 ARTICULATOR

Physical and virtual articulator designed with special geometries to simulate patient-specific rotating, sliding and closing jaw movements. In Zirkonzahn scanners, models can be scanned together with the PSI and then transferred into the virtual articulator.

GS1 PLASTER ARTICULATOR

Articulator for plastering models and checking opening and closing movements, specially conceived to protect the PSI articulator from plaster residues and water. A specific calibration key creates a synchronisation with the PSI.

If desired, this articulator can also be used in combination with the JawAligner system, without using plaster.

MINI-ARTI ZS1

Small articulator conceived for checking the occlusion of small bridges and single crowns, which can also be used for full-arch models. Functional and easy to use, it has been designed for prototyped models (milled or printed) and as valuable instrument for simulating opening, closing, protrusion, retrusion and laterotrusive movements.









DIFFERENT HEIGHTS

The JawAligner are available in different heights, permitting the articulation of models independently of their height. Depending on model requirements, the different JawAligner heights can be combined with each other to save resin during printing.



INCLINED AND STRAIGHT

The JawAligner system consists of two sets, one inclined and one straight, which allow to articulate models according to the patient-specific occlusal plane inclination. These two different sets were specially developed to provide models with great stability, avoiding deformations caused by an unbalanced front-rear height and permitting their articulation regardless of the inclination degree. Thanks to the different heights and inclinations of the two JawAligner sets, any model can be placed in the articulator, even, for example, implant models with particularly long laboratory analogs.



JawAligner PS1 Straight Low



JawAligner PSI Straight Medium



JawAligner PS1 Straight High



JawAligner PS1 Inclined Low




JawAligner PS1 Inclined High



MODELS CREATION WITH THE JAWALIGNER SYSTEM

Especially when dealing with printed models, the hardening process of plaster might cause deformations. The innovative JawAligner system, completely integrated into Zirkonzahn's digital workflow, not only ensures a plaster-free articulation but saves also time and material.



After taking a digital impression of the oral situation with the Detection Eye intraoral scanner, the acquired data can be transferred into the Zirkonzahn. Modifier software for models design.



To use the JawAligner, first of all choose the articulator in the Model Maker module of the Zirkonzahn. Modifier software.




Then choose the JawAligner set and the most suitable spacer plate, depending on model size and articulator used (PSI, or other articulators from other manufacturers).



Depending on model requirements, the spacer plates can be combined with each other to save resin during printing.



Once the creation of the models is complete, the software creates the corresponding attachments in the model bases.



You can also set parameters between spacer plates and models (e.g. creation of a honey comb structure for more stability, increase of base thickness, attachment adaptation or compensation of possible layer imprecisions).

Zirkonzahn



At this point, models can be printed with the 3D P4000 System or milled in Zirkonzahn's milling units.



Fix the threaded screw up to its halfway point.



Remove the screw and clean the screw housing from all resin shavings left from the screwing process.



Note: according to the model weight, two different metal plates can be used, showing lower or higher magnetic strength:

1 mm and 1.5 mm.



Place the magnetic plate in the attachment and fix the screw. If desired, it is also possible to glue the plate directly to the model.



If there is any resin excess in the two small antirotation attachments, use a scalpel to remove it.



The final models with magnetic plates can be now plaster-free articulated.



Models are finally attached to the JawAligner in the articulator to simulate patient's jaw movements.



IMPROVED COMMUNICATION

Dental models designed with the JawAligner system can be articulated in both PSI and Mini-Arti ZSI articulators, as the heights of the JawAligner ZSI recreates the same of the JawAligner PSI Inclined. The Mini-Arti ZSI is the ideal tool a dental technician can use to deliver or discuss with the dentist any kind of dental job produced from intraoral scanners. This means this articulator stands out not only for its suitability, but also for being a useful professional tool for improved clinic-laboratory communication.







CALIBRATION OF THE ARTICULATORS

A correct calibration ensures the right distance between the upper and lower split-cast, allowing the correct transfer of the digitally recorded occlusion into the laboratory articulator. For calibrating the articulators, specific kits are available, containing Adjusting Key and corresponding Zero Adjusting Pin.

HOW TO CALIBRATE THE PS1 ARTICULATOR

To use the JawAligner with the PSI articulator, this has to be previously calibrated as follows in order to achieve maximum precision:

- 1. Insert the Adjusting Key PSI and the Zero Adjusting Pin PSI into the PSI articulator without the upper split-cast.
- 2. Place the upper split-cast on the calibration key, apply a glue layer and close the articulator.
- 3. To avoid any distorsion while the glue is hardening, place a light weight on the anterior part of the articulator.
- 4. After the glue hardening is complete, unscrew the screw which joins the two parts of the calibration key and open and close the articulator, checking that the two parts of the Adjusting Key PSI fit together correctly.



INSTRUCTION: CALIBRATION OF THE PS1 ARTICULATOR



Zero Adjusting Pin PS1 + Adjusting Key PS1



COMPATIBLE WITH ARTICULATORS OF OTHER MANUFACTURES

WITH THE CORRESPONDING ACCESSORIES, THE JAWALIGNER ARE ALSO COMPATIBLE WITH ALL ARTEX ARTICULATOR OF 126 mm HEIGHT, WITH SAM 2P, 2PX AND 3 ARTICULATORS (NOT WITH +15 VERSIONS) AND KAVO.



ALL ARTEX ARTICULATORS OF 126 MM HEIGHT



All Artex articulators of 126 mm height



Set Adjusting Key Artex Item number: ZBAC3611

The range of compatible articulators is constantly expanded.

SAM 2P/2PX/3 ARTICULATORS (NOT WITH +15 VERSIONS)



SAM 2P/2PX/3



Set Adjusting Key SAM Item number: ZBAC3621



Set JawAligner PS1 Adapter SAM Item number: ZBAC3622

KAVO



KaVo Protar evo 7/evo 9



Set Adjusting Key KaVo Item number: ZBAC3631



Set JawAligner PS1 Adapter KaVo Item number: ZBAC3632



MODELS IN THE MINI-ARTI ZS1

If models have been designed using the PSI virtual articulator and the Set JawAligner PSI Inclined, they can also be fixed in the Mini-Arti ZSI once milled or printed by means of the Set JawAligner ZSI.





Set JawAligner ZS1: two magnetic spacer plates (High and Low) for each dental arch, screws, Allen key and magnetic plates



Step 1: screw the Upper Jaw Low spacer plate at the top of the Mini-Arti ZS1



Step 2: screw the Lower Jaw Low spacer plate at the bottom of the Mini-Arti ZSI



Step 3: if models have been designed in the software using the High spacer plates of the JawAligner PSI Inclined, attach in the Mini-Arti ZSI the respective additional plate (Upper Jaw High/Lower Jaw High)



Step 4: at this point, fix the models into the Mini-Arti ZS1 ...



Step 5: ... and check the restoration occlusion



OVERVIEW OF THE JAWALIGNER SETS



SET JAWALIGNER PS1 INCLINED

- 4 magnetic spacer plates (2 x dental arch) with an inclination of 5° and two differents heights (High and Low)
 - Upper Jaw High
 - Upper Jaw Low
 - Lower Jaw High
 - Lower Jaw Low
- Magnetic plates Ø 25 x 1,5 mm
- Magnetic plates Ø 25 x 1 mm
- Countersunk-head self-tapping screws TX10
- Allen key TX10
- Wooden toolbox

Item number: ZBAC3561



SET JAWALIGNER PS1 STRAIGHT

- 6 magnetic spacer plates (3 x dental arch) in three different heights (High, Medium and Low)
 - Upper Jaw High
 - Upper Jaw Medium
 - Upper Jaw Low
 - Lower Jaw High
 - Lower Jaw Medium
 - Lower Jaw Low
- Magnetic plates Ø 25 x 1,5 mm
- Magnetic plates Ø 25 x 1 mm
- Countersunk-head self-tapping screws TX10
- Allen key TX10
- Wooden toolbox

Item number: ZBAC3591



SET JAWALIGNER ZS1

- 4 magnetic spacer plates (2 x dental arch) for two different heights (High and Low)
 - Upper Jaw High
 - Upper Jaw Low
 - Lower Jaw High
 - Lower Jaw Low
- Magnetic plates Ø 25 x 1,5 mm
- Magnetic plates Ø 25 x 1 mm
- Countersunk-head self-tapping screws TX10
- Hexagon countersunk-head screws
- Hexagon Allen key
- Allen key TX10

Item number: ZBAC3571








ZIRKONZAHN WORLDWIDE

In our headquarters in Gais (South Tyrol, Italy) we take care of all tasks:
from research and development, to customer support for technical and dental
technical issues, all the way through the sale of our products. In order to offer
all our customers the best possible support as well as optimal access to our
training courses, we have established additional offices and education centres
in ten countries.



THE JAWALIGNER SYSTEM

Zirkonzahn Worldwide – An der Ahr 7 – 39030 Gais/South Tyrol T +39 0474 066 680 – F +39 0474 066 661 – www.zirkonzahn.com – info@zirkonzahn.com