











General instructions


Veneering of Tecno Med Mineral bridges with composites

Required tools and materials

- | | | | | | |
|--|--|---|-------------------------|---|--------------------------------------|
|  | Sand blaster (aluminium oxide: grain size 110 µm; 2,5 bar) |  | Light curing lamp |  | Fine, cross-cut, hard-metal burs |
|  | Dry and oil-free compressed air ISO/DIN 7494-2:2014-03 |  | Opaque* and brush |  | Polishing paste* and polishing brush |
|  | Ethanol* (without parabens and scent-free) and disposable applicator |  | Composite Veneer resin* | | |
|  | Composite Bonding Liquid for Tecno Med (CPAA2060)* and swab |  | Glycerine gel* and swab | | |


⚠ * ATTENTION! Please follow the respective instructions from the manufacturer! Differences are possible due to simplified descriptions.

Pre-treatment/cleaning of the Tecno Med Mineral bridge



- 1 Sand blast and clean with compressed air
- 2 Degrease with ethanol*
- 3 Apply the Composite Bonding Liquid for Tecno Med*
- 4 Polymerise (180 secs.)
- 5 Apply the opaque*
- 6 Polymerise (180 secs.)


Modelling of the bridge with composite



- 1 Apply composite* to form the tooth details (max. 2 mm)
- 2 Polymerise (180 secs.)
- 3 Apply the last layer of composite*
- 4 Coat with glycerine gel*
- 5 Polymerise (180 secs.)

If necessary, repeat the steps 1 – 2 until you get the desired result. Use different dentine coloured composites and polymerise them after each layer.

Finishing of the bridge



- 1 Elaborate with a fine, cross-cut, hard-metal bur
- 2 Buff the surface with polishing paste*
- 3 Finished bridge

⚠ * ATTENTION! Please follow the respective instructions from the manufacturer! Differences are possible due to simplified descriptions.