


# Zirkonzahn®

*Human Zirconium Technology*



**100% PRETTAU ZIRCONIA - MADE BY CAD/CAM**

*“Perfection through perseverance“ Enrico Steger*



*Zirkonzahn World Wide*



*Zirkonzahn Germany*



*Zirkonzahn USA*




*Zirkonzahn Mexiko*

## *Prettau Zirconia*

In order to demonstrate the functionality of these types of construction, we collaborate with Dr. Fernando Rojas-Vizcaya from Spain and manufacture several works each month. The advantages of occlusally screwed Prettau bridges are their easy fabrication and their stability.

Using wax-up and double-scan, the construction is realized in a short period of time. The system creates the occlusal screw accesses automatically. The buccal part of the front teeth has been veneered, resulting in a very high esthetic in the visible part and high stability in the lateral tooth area.



## *The Prettau Bridge*

Dr. Fernando Rojas-Vizcaya, Spain

Georg Walcher – Zirkonzahn Education Center




# Zirkonzahn®



## HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN World Wide - Tel +39 0474 066 680 - [info@zirkonzahn.com](mailto:info@zirkonzahn.com) - [www.zirkonzahn.com](http://www.zirkonzahn.com)



WEAA5429=

*Milled with Zirkonzahn CAD/CAM System 5-TEC*