

Zirkonzahn®

Human Zirconium Technology




100% PRETTAU ZIRCONIA

„Perfection through perseverance“ Enrico Steger

Prettau Zirconia

We decided to restore this full upper arch with a 'cement on' full-zirconia Prettau bridge. All abutments were milled individually, two-tone coloured and fitted with titanium connectors to provide extra stability. The lower front segment was retained by lingual screws. With the best aesthetic outcome in mind vertical bone loss was compensated for by tissue coloured flanges veneered with pink ICE Tissue Ceramic.

The anatomical tooth shapes were milled in Prettau zirconia entirely (no veneer porcelain added!)



The Prettau Bridge

Clinician: Dr. Ghadeer Thalji – University of North Carolina, Department of Prosthodontics

Technician: Georg Walcher – Zirkonzahn Education Center




Zirkonzahn®



HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN World Wide - Tel +39 0474 066 680 - info@zirkonzahn.com - www.zirkonzahn.com



Milling without limits