



Zirkonzahn®

Human Zirconium Technology



BITE SPLINTS

Software-module CAD/CAM











SOFTWARE-MODULE CAD/CAM BITE SPLINTS

Module for the manufacture of bite splints

- *precise modelation of the necessary contact points in conjunction with the virtual articulator*
- *individual shaping*
- *adjustable parameters: occlusal thickness, wall thickness, smoothness (blocking-out)*
- *millable in various plastic materials*

Item number: SY0530





Zirkonzahn®



HUMAN ZIRCONIUM TECHNOLOGY

Zirkonzahn Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com

Milled with the M1 Wet Heavy Metal CAD/CAM system



WEAA5493=

All information is subject to change. Errors and omissions excepted. Version: 20/03/2018