

**Zirkonzahn<sup>®</sup>**

*Human Zirconium Technology*



**PRETTAU<sup>®</sup> ANTERIOR<sup>®</sup>**



Zirkonzahn®



*Full contour structures made from Prettau® Anterior® and colored with ICE Zirkon Malfarben 3D by Enrico Steger Stains*



# PRETTAU® ANTERIOR®

*Prettau® Anterior® is a new and highly biocompatible zirconia material. Due to its excellent esthetical properties and higher transparency, it is ideally suited for use in the anterior tooth region. Moreover, it is also perfectly appropriate for use in the posterior tooth region.*

- *Higher transparency due to optimized microstructure of the zirconia*
- *Flexural strength: 670 MPa (School of Dental Medicine, University of Pennsylvania)*
- *Indications: metal-free partial and single crowns, inlays, onlays, veneers, max. 3-unit bridges*
- *Fully anatomical or reduced structures for ceramic layering*
- *Shorter sintering time due to final sintering temperature of 1,500 °C (2,732 °F)*
- *For dry processing with common zirconia burs*
- *Customizable with Colour Liquid Prettau® Anterior® Aquarell, ICE Zirkon Keramik Ceramics, ICE Zirkon Malfarben Prettau® Stains and ICE Zirkon Malfarben 3D by Enrico Steger Stains*



# Zirkonzahn®



*The perfect anterior tooth esthetics with Prettau® Anterior® colored with ICE Zirkon Malfarben 3D by Enrico Steger Stains*

## HUMAN ZIRCONIUM TECHNOLOGY

*ZIRKONZAHN Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com*

*ZIRKONZAHN USA Inc. – Phone +1 800 989 8931 – Phone +1 678 441 9419 – infousa@zirkonzahn.com*



*Prettau® Anterior®*



*Prettau® Zirconia*