

Zirkonzahn<sup>®</sup>

*Human Zirconium Technology*



**PRETTAU<sup>®</sup> ANTERIOR<sup>®</sup>**



## INCREASE YOUR PRICE WITH PRETTAU® ANTERIOR®

**Transparency combined with high flexural strength!** The new zirconia material has been specifically developed for anterior tooth restorations, but is just as well suited for posterior areas. Prettau® Anterior® can be milled dry with regular zirconia burs and therefore offers an ideal alternative to lithium disilicate!

- Double the flexural strength compared to lithium disilicate: 670 MPa (School of Dental Medicine, University of Pennsylvania)
- Excellent esthetic properties
- Higher transparency allows perfect light transmission and good background masking
- Extra resistant to ageing
- Suitable for fully anatomical contours or framework for ceramic veneering
- Shorter sinter times due to sintering at 1,500 °C (2,732 °F)
- Available in various block sizes
- Indications: metal-free partial and individual crowns, inlays, onlays, veneers, 3-unit bridges (max).





# Zirkonzahn®



*Perfect esthetics also in the posterior area with Prettau® Anterior®, colored with ICE Zirkon Malfarben 3D by Enrico Steger surface stains.*

## HUMAN ZIRCONIUM TECHNOLOGY

*ZIRKONZAHN Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com*

*ZIRKONZAHN USA Inc. – Phone +1 800 989 8931 – Phone +1 678 441 9419 – infousa@zirkonzahn.com*



*Prettau® Anterior®*



*Prettau® Zirconia*